

Exhibit I

YAHOO! FINANCE [Home](#) - [Yahoo!](#) - [Help](#)**PR Newswire**


[[Business](#) | [US Market](#) | [By Industry](#) | [IPO](#) | [AP](#) | [S&P](#) | [International](#) | [PRNews](#) | [BizWire](#) | [CCN](#)]

Monday June 7, 7:57 am Eastern Time

Company Press Release

SOURCE: SMLX Technologies, Inc.

SMLX Technologies, Inc.'s Joint Venture With HelveStar, S.A., Creates Global Force For Core Product Distribution



HALLANDALE, Fla., June 7 /PRNewswire/ — SMLX Technologies, Inc. (OTC Bulletin Board: **SMLX** - **news**) is building a global force for distribution of its key products. The company's recent joint venture with HelveStar, S.A., a Swiss technology, commercialization, financial and holding company, formed BioStar Life Sciences, which will have a 99-year worldwide exclusive license to produce, market and distribute all SMLX core products and technology. The deal gives SMLX a 40-percent stake in the Switzerland-based BioStar, and the right to name two members to BioStar's board of directors.

BioStar's opening order included 300,000 Airbrators, 42,500 doses of Bioven and 300,000 Rapid Saliva Test kits. BioStar also purchased certain non core products from SMLX for \$2.4 million, payable \$100,000 per month for 24 months, and will underwrite the cost to commercialize these products.

HelveStar, engaged in pharmaceutical, nutraceutical, cosmeceutical and medical diagnostics, will own 60 percent of BioStar and have the right to name two directors to SMLX's board of directors. HelveStar, during the next 24 months, can purchase up to 2,250,000 shares of SMLX stock and 900,000 shares through stock options.

* HelveStar has already taken key steps to implement BioStar's global distribution strategy. HelveStar has purchased SMLX's offshore contract manufacturer, and also bought control of a U.S. dental distributor to sell the Airbrator in the U.S. The Malaysian company, Usra Sdn Bhd, will begin manufacturing the Airbrator in Asia early next year. AsiaStar Healthcare will distribute the product in Asia.

BioStar is making plans to launch the sale of Bioven, SMLX's treatment for Rheumatic Arthritis, in Hong Kong, Malaysia and Singapore starting in July 1999. Also in July 1999, BioStar plans to begin selling SMLX's diagnostic kits in several countries. In June, several key personnel from HelveStar's Malaysian companies will visit SMLX to finalize further details of the transfer of responsibility.

SMLX is a bio-engineering company that develops leading-edge, high-technology products for the Life

BioStar Airbrator

— FOR PROFESSIONAL USE ONLY —

Order Form

Name/Company _____

DDS / DMD

Address _____

City _____

State _____

Zip _____

Phone () _____

Ordered by _____

Specialty _____

Item	Qty.	Price Each	Qty.	Price Per Box of 10	Amount
<input type="checkbox"/> High Performance Airbrators*	_____	\$14.95	_____	\$99.95	_____
<input type="checkbox"/> Medium Performance Airbrators*	_____	\$14.95	_____	\$99.95	_____
<input type="checkbox"/> Light Performance Airbrators*	_____	\$14.95	_____	\$99.95	_____
<input type="checkbox"/> Airbrator* Starter Kit	_____	\$84.00	_____	_____	_____
<input type="checkbox"/> Connector	_____	\$35.00	_____	_____	_____
<input type="checkbox"/> Tip Bending Tool	_____	\$5.00	_____	_____	_____

☐ Bill Me

☐ Payment Enclosed

Make checks payable to:
Edge Dental, Inc.

Sub-Total _____

MI Residents add 6% Sales Tax _____

Shipping & Handling _____


Total _____

Call _____

Or order by phone:

1-800-873-6070

Fax 231-922-2274




CLIP & SEND

Call today for a
free catalog of products from
Edge Dental including:

- Diack/Vac — Chemical Indicators
- Sterilization Pouches
- Palatal Anesthesia Device
- Edgemate Sharpening Stones
- Disposable Chairside Sharpening Guide
- Hartzell Dental Instruments
- NorthBay/Bioscience Sterilizer Monitoring Service


1-800-873-6070



AIRBRATORS® \$14.95 each
Box of 10 Airbrators®: \$99.95 per box
Tip Bending Tool: \$5.00
Connector: \$35.00

Pricing

STARTER KIT \$84.00 includes:
3 High Performance Airbrators®
1 Medium Performance Airbrator®
2 Light Performance Airbrators®
1 Airbrator® Connector, and Instructions



Bulk Rate
U.S. Postage
PAID
Grand Rapids, MI
Permit #657

EDGE DENTAL, INC.
PO BOX 4332
TRAVERSE CITY, MI 49685-4332


www.edgedental.com

Call for a free catalog
1-800-873-6070

EDGE

BioStar Airbrator

- Patented delivery system maintains constant particle flow, increasing operator control and minimizing powder overspray.
- Disposable handpiece reduces the threat of cross-contamination.
- Quick connect to existing dental unit high speed air line minimizes setup and cleanup time for dental assistants.
- Handpiece design requires no maintenance, no sterilization, no filling with messy powders, no expensive tip replacements, and doesn't require additional floor or countertop space.
- FDA-registered.



Connector screws onto existing air line.

What are the advantages of the Airbrator over traditional air abrasion units?

The major advantage is cost. Traditional air abrasion units can cost up to \$20,000! Airbrators* cost as little as \$9.95 each, with these added advantages:

- No maintenance, no refilling, no autoclaving;
- No dragging heavy equipment from operator to operator;
- Fits standard handpiece connectors; and
- Setup is simple.

Will I be able to use my new Airbrators immediately?

Yes, however, as with other new dental equipment, some practice and study of technique is strongly recommended. We recommend particular emphasis on case selection.

For what procedures should I use the Airbrator?

Applications include class I, III, IV & V preparations, removing composites, porcelain and amalgam stains and P.R.R. (Preventive Resin Restoration) preparations.

Does the Airbrator require any special care?

No. Simply remove from its pouch, connect to your air line, use, and dispose. No need to autoclave. No need to refill.



Use of a rubber dam is recommended.

Will the Airbrator save me time?

The major benefit of air abrasion is the lack of sensitivity to the patient and less discomfort for the patient. The reduced need for anesthesia, multiple quadrant procedures and minimal set up time all add up to a more productive, time

Air Abrasion Technology

- Conserves the maximum amount of healthy tooth structure.
- Saves time. Doctors report that less than 10% of their patients need anesthetic when air abrasion is used, saving injection and numbing time.
- Allows for multiple quadrant treatment in one visit.
- Quiet. No vibration, and no high-pitched noise.
- Less sensitivity to your patients.
- Patients love it! Procedures take less time. Patients experience less postoperative sensitivity.



Tip bender reduces crimping.

How can I minimize overspray?

As with traditional air abrasion units, for intraoral applications you should use high speed suction. For patient protection, a rubber dam or moistened gauze pads surrounding the work area and protective eyewear should be used. For extra-oral applications, use a dust cabinet. Although overspray is a common complaint with all types of air abrasion, practitioners commonly state that overspray can be greatly diminished with experience.

I'm not sold yet. Why else should I give the Airbrator a try?

Wouldn't your patients like to know that you take a less invasive, more preventive approach to dentistry, treating early stage caries rapidly and comfortably rather than watching and waiting? Don't

BioStar® Airbrator

- Patented delivery system maintains constant particle flow, increasing operator control and minimizing powder overspray.
- Disposable handpiece reduces the threat of cross-contamination.
- Quick connect to existing dental unit high speed air line minimizes setup and cleanup time for dental assistants.
- Handpiece design requires no maintenance, no sterilization, no filling with messy powders, no expensive tip replacements, and doesn't require additional floor or countertop space.
- FDA-registered.



Connector screws onto existing air line.

What are the advantages of the Airbrator over traditional air abrasion units?

The major advantage is cost. Traditional air abrasion units can cost up to \$20,000! Airbrators® cost as little as \$9.95 each, with these added advantages:

- No maintenance, no refilling, no autoclaving;
- No dragging heavy equipment from operator to operator;
- Fits standard handpiece connectors; and
- Setup is simple.

Will I be able to use my new Airbrators immediately?

Yes, however, as with other new dental equipment, some practice and study of technique is strongly recommended. We recommend particular emphasis on proper technique.

For what procedures should I use the Airbrator?

Applications include class I, III, IV & V preparations, removing composites, porcelain and amalgam stains and P.R.R. (Preventive Resin Restoration) preparations.

Does the Airbrator require any special care?

No. Simply remove from its pouch, connect to your air line, use, and dispose. No need to autoclave. No need to refill.




Use of a rubber dam is recommended.

Will the Airbrator save me time?

The major benefit of air abrasion is the lack of sensitivity to the patient and less discomfort for the patient. The reduced need for anesthesia, multiple quadrant procedures and minimal set up time all add up to a more productive, time

Air Abrasion Technology

- Conserves the maximum amount of healthy tooth structure.
- Saves time. Doctors report that less than 10% of their patients need anesthetic when air abrasion is used, saving injection and numbing time.
- Allows for multiple quadrant treatment in one visit.
- Quiet. No vibration, and no high-pitched noise.
- Less sensitivity to your patients.
- Patients love it! Procedures take less time. Patients experience less postoperative sensitivity.



Tip bender reduces crimping.

How can I minimize overspray?

As with traditional air abrasion units, for intraoral applications you should use high speed suction. For patient protection, a rubber dam or moistened gauze pads surrounding the work area and protective eyewear should be used. For extra-oral applications, use a dust cabinet. Although overspray is a common complaint with all types of air abrasion, practitioners commonly state that overspray can be greatly diminished with experience.

I'm not sold yet. Why else should I give the Airbrator a try?

Wouldn't your patients like to know that you take a less invasive, more preventive approach to dentistry, treating early stage caries rapidly and comfortably rather than watching and waiting? Build

BioStar™ Airbrator®

POLISHING & CLEANING

See page insert for instructions

Call Now! 1-800-873-6070

Caution: This is a disposable device. Please do not reuse. Do not attempt to sterilize, re-sterilize, or reuse. For more information, please call 1-800-873-6070 or visit our website at www.edge-dental.com.



Manufactured under exclusive license for BioStar Life Science, SA by:
EDGE DENTAL, INC.
 13606 S. West Bay Shore Drive
 Traverse City, MI 49684
 Telephone 1-800-873-6070
 Fax (231) 928-2274

Simplex® Airbrator®

HIGH PERFORMANCE ABRASION

Manufactured by SMLX Technologies, Inc.

378 Arsin Blvd., Hallandale, FL 33009
 Telephone (954) 455-0110
 Fax (954) 455-9008



High Performance Abrasion
 Aluminum Oxide - 5.0g
 For Professional Use Only

Medium Performance Abrasion

Aluminum Oxide — 5.0g

For Professional Use Only

BioStar™ Airbrator®

MEDIUM PERFORMANCE ABRASION

See package insert for instructions.

Cat. No.: 100-302 Lot No.: 1000

Caution: This is a disposable device. Discard after single patient use. Do not attempt to sterilize, refill or reuse. Federal law restricts this device to sale by or on the order of a dentist.

Manufactured under exclusive license for BioStar Life Sciences, SA by:
SMLX Technologies, Inc.
 376 Ansin Blvd., Hallandale, FL 33009
 Telephone (954) 455-0110
 Fax (954) 455-9008

Distributed by:
EDGE
 Edges Dental, Inc.
 13606 S. West Bay Shore Drive
 Traverse City, MI 49684
 Telephone 1-800-873-6070
 Fax (231) 928-2274

