


Exhibit D

Assembly Step #3: The Inlet Tube is press fit into the Body from the Top a specified depth.



Assembly Step #5: Inlet Seal is pressed into the Body from the Top a specified depth, so it makes contact with the bonding agent. Bonding agent wells around both parts to the body and then cured.

Assembly Step #4: The bonding agent is dispensed along the peripheral of the Inlet Tube from the Top of the Body.



KUS Products, Inc.
954-476-7736
954-764-4074
Dental Convention for Dispensing
Full Seal
Code 471069